



# Instructions to install Python and Anaconda on a Windows-running machine.

The instructions assume having access to a Microsoft Visual Studio installer.

1. Execute Microsoft Visual Studio installer and click on individual components.




2. Select and install Python 3 and Anaconda (32 or 64 bit, depending on the machine).




3. PyDicom Library installation.

For this step it is necessary to run Anaconda as administrator.

First, locate Anaconda in the startup menu and open the file location (right click on the shortcut):



Then, run Anaconda as administrator (right click on the shortcut):



In the Anaconda Prompt terminal, run “pip install pydicom”:

```
(base) C:\Windows\system32>pip install pydicom
Collecting pydicom
  Downloading https://files.pythonhosted.org/packages/43/88/d3c419ab2e753e7651510882a53219373e78fb55294cb247dffd3934ea55/pydicom-1.2.2-py3-none-any.whl (7.0MB)
100% |#####| 7.0MB 2.0MB/s
Installing collected packages: pydicom
Successfully installed pydicom-1.2.2
You are using pip version 10.0.1, however version 19.0.2 is available.
You should consider upgrading via the 'python -m pip install --upgrade pip' command.
(base) C:\Windows\system32>
```

Once the process is completed, the PyDicom library is installed on the system and ready to use.

## Instructions to run .py scripts

- After downloading or cloning the TeleMed repository, open up Anaconda (no Administrator privilege needed), navigate to the project’s folder, and run the script/application following the instructions inside of the script.

